

Blood Pressure Monitor Performance Report

Customer:	KU-CTSU Rainbow	Contact:	Andra Lahner
Street:	3901 Rainbow Blvd (Sudler 2nd Floor)	Phone #	913-574-3018
City/ST/Zip:	Kansas City, KS 66160	P= Pass F=Fail	

MFG	Model#	S/N	Electrical Safety Check	Condition Controls/ Switches	Condition Power Cord & Strin Relief	Hardware Integrity	Indecators & Display Operation	MFG. Labeling Intact	Normal		Hypertonic		Hypotonic		SPO2	Date Performed
									120/80	P 60	200/150	P130	60/30	P40	97%	
AXIA Surgical	TRIA	A18014600076	P	P	P	P	P	P	120/83	59	203/154	130	60/31	39	97%	1/19/2021
AXIA Surgical	TRIA	A18044501442	P	P	P	P	P	P	122/83	60	203/152	130	62/32	40	98%	1/19/2021
AXIA Surgical	TRIA	A18044501468	P	P	P	P	P	P	120/83	60	205/151	130	60/31	40	96%	1/19/2021
DELPG039																
Welch Allyn	297	200819211	P	P	P	P	P	P	119/82	60	201/152	131	60/33	40	NA	1/19/2021
Welch Allyn	42NTB	201907174	P	P	P	P	P	P	121/81	60	202/153	130	60/32	40	98%	1/19/2021

Technicians	JC	Technician:	Jim Cowart
Initials		Date:	01/19/2021
		Test equipment used	
		FLUKE- ProSim 8 S#399003 Expires:09/11/2021	

Notes: Next Due Date: 01/2022

Calibration Performed by Cobb Refrigeration and Lab Service 8071 Flint St. Lenexa KS, 66212 913-403-0088

Blood Pressure Monitor Performance Report

Customer: KU-CTSU Rainbow		Contact: Andra Lahner														
Street: 3901 Rainbow Blvd (Sudler 2nd Floor)		Phone #: 913-574-3018														
City/ST/Zip: Kansas City, KS 66160		P= Pass F=Fail														
MFG	Model#	S/N	Electrical Safety Check	Condition Controls/ Switches	Condition Power Cord & Strin Relief	Hardware Integrity	Indecators & Display Operation	MFG. Labeling Intact	Normal		Hypertonic		Hypotonic		SPO2	Date Performed
									120/80	P 60	200/150	P130	60/30	P40		

Blood Pressure Monitor Performance Report

Customer: KU-CTSU Rainbow		Contact: Andra Lahner														
Street: 3901 Rainbow Blvd (Sudler 2nd Floor)		Phone #: 913-574-3018														
City/ST/Zip: Kansas City, KS 66160		P= Pass F=Fail														
MFG	Model#	S/N	Electrical Safety Check	Condition Controls/ Switches	Condition Power Cord & Strin Relief	Hardware Integrity	Indecators & Display Operation	MFG. Labeling Intact	Normal		Hypertonic		Hypotonic		SPO2	Date Performed
									120/80	P 60	200/150	P130	60/30	P40		